




DSC Testing Report for Ethylene Oxide

Report #:	ECUST-RCSC-20240419-O01		
Institute:	Reactive Chemical Safety Center at East China University of Science and Technology, Shanghai, China. RCSC@ECUST (https://rcsc.ecust.edu.cn/)		
Analyst:	Shiyi Li	Test time:	2024-04-19
Reporter:	Shiyi Li	Report time:	2025-07-15
Reviewer:	Min Sheng	Review time:	2025-07-15
Limits of Applicability:	The testing data generated by the RCSC@ECUST is based on specific samples from a certain supplier. During the testing process, we have made every effort to ensure the reliability of the data, but we still cannot exclude data quality issues caused by supplier differences, sample impurities, deterioration of long-term storage, testing method error, and human error. Therefore, we suggest only using RCSC testing data as a reference, and we do not assume any responsibility for any losses caused by using this data.		
Data Copyright:	Prohibit the sharing of this report on websites outside the Reactive Chemical Safety Center at East China University of Science and Technology. If you have any questions, please contact RCSC@ECUST (https://rcsc.ecust.edu.cn/).		

Chemical Names	Ethylene Oxide, EO
Synonyms	Epoxyethane 1,2-Epoxyethane Oxacyclopropane Dihydrooxirene Oxidoethane
Molecular Formula	C2H4O
Hill Notation Formula	C2H4O
CAS #	75-21-8
Test	DSC_N2
Notes	Glass capillary, 0~500 °C, 0.8°C/min, Nitrogen headspace



Detail information of the tested sample:

Project	Development of Reactive Chemical Database for RCSC@ECUST
Sample Description	
Composition	99.5 %
Chemical Structure	
Supplier	Alladin
Storage Information	
Sample Status	Transparent liquid
Other notes	

Testing Result:

1 st test					
Instrument: TA DSC Q2500		DSC crucible: Glass Capillary		Headspace gas: N ₂	
Tested sample mass: 1.114 mg		Total pre-test mass: 28.216 mg		Total post-test mass: 28.216 mg	
Start Temperature: 0°C		Scan rate: 0.8°C/min (Twice)		End Temperature: 500°C	
#	Endo/Exo Peak	Onset T (°C)	Peak T (°C)	End T (°C)	Total Heat (J/g)
1	Exotherm	308	427	472	-2792.8 ^a
Notes: ^a , TNT equivalent 0.67					




Figure 1. DSC Curve with Auto scale, the solid line represents the first scan (with integration); The dashed line represents the second scan




Figure 2. DSC Curve with -6~6 W/g scale, the solid line represents the first scan (with integration); The dashed line represents the second scan



Figure 3. Tested Sample Pictures (original package, tested sample before test and after test)